

sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

edition 2024/01

TIMS

BROCHURE

& ORDERING GUIDE

TECNOLOGIA INOVADORA MAIS SEGURA, LDA

TECHNOLOGICAL INNOVATION WITH SAFETY, LDA

TIMS

Technological Innovation With Safety

© 2024 TIMS. Ed. 2024/01

TIMS. Ordering Guide

TIMS streamlined the ordering process into three simple steps to ensure you a smooth and efficient experience

Review our catalogue

Begin by reviewing our comprehensive catalogue, *which includes detailed information about each of TIMS's medical devices.*
In the catalogue, you will find product specifications, advantages, and dimensions, providing everything you need to make an informed decision.

Gather information

After selecting the desired products, **take note of the product references and intended quantities**. *This information is crucial for placing an accurate order and ensuring we deliver exactly what you need.*

Place your order

To finalize your order, scan the QR code provided, which will lead you to an **email form**. *Input the gathered product information, along with invoicing details and the shipping address. If you are an existing client, please include your institution's details for accurate processing. Alternatively, you can also send your order directly* **info@tims.pt** *email address for convenience.*

SCAN ME

We kindly ask that you **include a direct telephone/ mobile contact amd your address in your email,** *should we need to clarify any details or answer further questions regarding your order.*
We look forward to assisting you with your medical device needs.

IMAGE DISCLAIMER

- Images shown in this brochure may differ slightly from the actual products.
- Illustrations are used in place of real-life images when they provide a clearer or more effective representation of the product's use.

SIZE CHART

SYMBOLS & MARKS

Product-specific measurements are provided whenever possible. For convenience, a general size chart has also been developed to assist in selecting the appropriate size.

S	Small	M	Medium	L	Large	XL	Extra Large
XH	Extra Height	LP	Low Profile				

ACRONYMS

For ease of reading and clarity, certain acronyms are utilized throughout the text to streamline the content.

MF - Memory Foam

TIMS PRODUCT SPREAD

NEW PRODUCTS & UPDATES

TIMS

Edition 2024/01 © 2024 TIMS. This brochure is regularly reviewed and updated to reflect our expanding product range. As new innovations and products are introduced, the brochure will be revised to ensure you have access to the latest information and offerings. Please refer to the most current version for accurate and up-to-date details.
If you need a device that is not currently available, please talk to us. Our designers may be able to collaborate with you and the clinicians to create a new product.

COPYRIGHT

© 2024 TIMS. All rights reserved. The content, images, and illustrations contained in this brochure are the exclusive property of TIMS and have been either captured or commissioned by the company. Unauthorized reproduction, distribution, or use of any material is strictly prohibited without prior written consent from TIMS.

INDEX.

TIMS BROCHURE EDITION DETAILS

Image Disclaimer 4
Symbols & Marks 4
TIMS Product Spread: New products & updates 4
Copyright 4

MEMORY FOAM

TIMS Memory Foam Product Selection 6

SAFER PRONE POSITIONING

Prone Head Supports & Memory Foam Cushions 8
Torso Support 22
Dome Supports: Chest & Pelvic 28
Arm Supports 38
Leg Supports 41

SAFER SUPINE POSITIONING

Supine Head Support / Beach Positioning 47

SAFER PATIENT CARE

Safer Patient Positioning
- Protective Straps 50
- Mirrors 52
- Helmet Angle Pads 52

Safer Eye Protection
- Safer Clinician Eye Protection 53
- Safer Patient Eye Protection 53

Safer Patient Health Protection
- Protective Mitts 54

© 2024 TIMS. Ed. 2024/01

sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

TIMS

SAFER SURGICAL POSITIONING

TIMS MEMORY FOAM PRODUCT SELECTION

All our toppers are supplied in memory foam. The viscoelastic properties of memory foam allow TIMS' Patient Support Cushions to mold to the individual shape of each patient, using the heat and pressure from the patient's body. This personalized contouring to each patient's size and weight distributes their weight more evenly, relieving pressure points and providing extra comfort without resistance. By evenly distributing pressure across the patient's body, muscles are able to relax more, preventing the build-up of tension in isolated areas. TIMS' design, incorporating memory foam's viscoelastic properties, helps prevent tissue markings and injuries, and significantly reduces edema.

Memory viscoelastic foam is also antibacterial and hypoallergenic, providing maximum comfort for patients and making it ideal for asthma sufferers, who make up 1 in every 12 adults.









sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Prone Positioning.

Prone Head Supports & Cushions

SAFER PRONE POSITIONING




PRONE HEAD SUPPORTS & MEMORY FOAM FACE CUSHIONS


TIMS MF Prone Face Cushions help the clinician achieve a Natural, Neutral Neckline for the patient for safer positioning. Single Use MF Face Cushion to be used with the TIMS Airway, Safe Pro Airway & Safe Lite Kits, as well as TIMS Face Support Base.

- Protect your patient from pressure related injuries to their face, eyes, forehead, cheeks and chin.
- Provide better comfort for your patient in the prone position.
- Different Shapes, Large & Small, supplied with Head Measuring Band and Secretion Pad.

Safer Prone Positioning. Prone Head Supports & MF Cushions




PRONE HEAD SUPPORT. SPECIFICATIONS AND AVAILABILITY SHOWN IN THE FOLLOWING PAGES.



MF FACE CUSHION

The depth and density of the disposable **TIMS Memory Foam Face Cushion** *means there is a noticeable reduction in facial oedema, substantial reduction in periorbital oedema and a lower risk of post-operative visual loss or stroke.*


FC-1S / 1L; FC-2S / FC-2L



EASY TURNING

Use the latex free **head strap** *as an 'extra hand' to help you support your patient's head during turning. The Prone Smart Head Strap should be used as an aid to your usual turning procedure and not as a replacement. Use the* **hand grip** *to ensure you have a firm hold on your cradle.*


STHS-01



SAFE & SECURE TUBING

Use **ET Tube Clip** *to keep tubing in place and stop it tangling.*

STHS-01



NATURAL NEUTRAL NECKLINE

Use our **custom adjustable cogs** *to adjust the height on the mirror platform.* **Use the cog's guide for greater accuracy.** *May be used without the mirror base platform for reduced height. The cradle can go as low as 16 mm, as long as it does not compromise the patient's neutral neckline.*

COGS-03

Safer Prone Positioning. Prone Head Supports & MF Cushions


TIMS Safe Head Support System and Reusable Accessories

- Have better control when turning a patient into the prone position by using the **unique hand grip** on the front of the cradle.
- Always have a level surface with the **adjustable 3 prong Mirror Base** to help achieve a neutral neckline.
- Each Cradle is supplied with **both sizes** of Memory Foam Prone Face Cushion.
- Use the **adjustable mirror** to monitor your patients eyes, endotrachael and nasal tubes. Prevent the patients tubing from becoming tangled by always using the ET Clip provided.



Easy cleaning for improved infection control. Clean TIMS' Safe Turn cradles with standard hospital wipes. The mirror can be removed for thorough cleaning.


(!) DO NOT AUTOCLAVE. DO NOT SUBMERGE IN WATER.



HEAD MEASURING BAND

TIMS' Head Measuring Band is designed to quickly and accurately determine the appropriate size of Prone MF Face Cushion for your patient – either Large or Small. The band features a straightforward color indicator that guides the physician in selecting the correct cushion to properly support the patient's face.

HMB-120-01




SECRETION PAD

Each kit comes with a Secretion Pad to keep your mirror base dry and free from liquids, helping to enhance infection control.

SEC-PAD-01

Safer Prone Positioning. Prone Head Supports & MF Cushions

SAFE TURN SYSTEM ACCESSORIES



FLAT & ANGLE PADS

Use TIMS Angle Pads to add extra height and a greater angle, if required.


HAP-01AD (ANGLE PAD SINGLE USE)
HAP-01-FLD (FLAT ANGLE PAD SINGLE USE)



MIRRORS

Monitor your patient's eyes, endotrachael and nasal tubes in the adjustable mirror platform. Use the small mirror to keep track of your patient's eyes when using the cradle directly on the table.


M-210X130



HEAD STRAP

Use the head strap to help you support your patient's head when turning them from supine to the prone position. The head strap should be used as an aid to your usual turning procedure, not as a substitute.


STHS-01



ETC CLIPS

Use the tubing clips to keep the tubing in place and stop it tangling.

ST-ETC-02



MIRROR BENEFITS


The mirror allows medical staff to ensure the patient's face remains correctly positioned during the procedure. This helps ensure that if there is any movement due to robust surgery, the patient's face can be safely repositioned, reducing the risk of pressure injuries. We recommend checking the patient's position every 30 minutes and documenting it.

TS-PRO-BASE-MC-01

Safer Prone Positioning. Prone Head Supports & MF Cushions


PRONE HEAD SUPPORT SYSTEMS

SAFE TURN PRO & AIRWAY



TIMS SAFE PRO

ST-PRO-CRADLE-01




TIMS AIRWAY

ST-AIRWAY-CRADLE-01

The TIMS Airway offers all the same features of the TIMS Safe Pro with one difference. The extra cut in the cradle allows you to position a patient without the need to disconnect the patient's ET Tubing (when using the MF Split Prone Face Cushion). Ideal for use in ICU and HDU when the patient is placed in the prone position.

Both kits include a mirror base with adjustable cogs.




MIRROR BASE W/ COGS

TS-PRO-BASE-MC-01

Safer Prone Positioning. Prone Head Supports & MF Cushions

PRONE HEAD SUPPORT SYSTEMS

SAFE LITE



TIMS SAFE LITE

ST-LITE-01-KIT

Safe Lite is our more affordable and portable system, equipped with a standalone mirror base to attach your head support cradle for surgeries performed with the patient in a prone position.

TIMS PRO, AIRWAY AND LITE HEAD SUPPORT SYSTEMS



ST-PRO-CRADLE-01 / ST-AIRWAY-CRADLE-01



ST-PRO-CRADLE-01 AND MIRROR BASE

Safer Prone Positioning. Prone Head Supports & MF Cushions


PRONE HEAD SUPPORT SYSTEMS

SAFE TURN KITS

references / SAFE TURN KITS AND COMPONENTS

TIMS Airway
ST-AIRWAY-01
2 per box


1 (pc) ST Cradle (Airway)
1 (pc) ST Mirror Base W/Cogs
1 (pc) ST Head Strap
2 (pcs) ST Loop Straps



ST-AIRWAY-CRADLE-01
ST-PRO-BASE-MC-01
STHS-01
STLS-01

TIMS Safe Pro
ST-PRO-01
2 per box


1 (pc) ST Cradle (Pro)
1 (pc) ST Mirror Base W/Cogs
1 (pc) ST Head Strap
2 (pcs) ST Loop Straps



ST-PRO-CRADLE-01
ST-PRO-BASE-MC-01
STHS-01
STLS-01

TIMS Safe Lite
ST-LITE-01
2 per box

1 (pc) Safe Turn Cradle (Lite)
1 (pc) Safe Turn Mirror Base
1 (pc) Mirror



ST-LITE-CRADLE-01
ST-LITE-BASE-01
M-210X130

All Kits include:

TIMS Safe Head Measuring Band	*THB-01*
TIMS Safe Storage Bag	*TSSB-01*
TIMS Safe Base (no mirror)	*TSB-01*
MF Face Cushion Large	*FC-1L (*)*
MF Face Cushion Small	*FC-1S (*)*
TIMS Secretion Pad	*SEC-PAD-01*
(*) Split Face Cushion alternative available	*FC-120S-SPLIT*
	FC-120-L-SPLIT

Safer Prone Positioning. Prone Head Supports & MF Cushions

PRONE HEAD SUPPORT SYSTEMS

MF PRONE FACE CUSHIONS



MEMORY FOAM PRONE FACE CUSHION FC 1S/1L; FC-2S/2L

Memory Foam Prone Face Cushion



MF PRONE FACE CUSHION DETAIL

Memory foam's viscoelastic properties allow the **TIMS MF Prone Face Cushion** to mold to the individual shape of each patient's face using the heat and pressure from the patient. This **personalized contouring** to the patient's face size and weight evenly distributes head weight, relieving pressure points and providing extra comfort without resistance. By spreading the pressure evenly across the forehead and cheek areas, facial muscles can relax more, preventing the build-up of tension in isolated areas. TIMS' design helps prevent tissue markings on the chin and reduces strain on the temporomandibular joints.

Memory foam is soft and provides more even support, whereas PU foam, traditionally used in prone face cushions over the last 25 years, can be abrasive on the patient's skin, particularly in older patients, younger patients with thin skin, and those on steroids.

Safer Prone Positioning. Prone Head Supports & MF Cushions


PRONE HEAD SUPPORT SYSTEMS

MF PRONE FACE CUSHIONS

Memory Foam Prone Face Cushion

The TIMS Prone Face Support Cushion's tilt design, combined with memory foam, supports correct spinal alignment and maintains the neutral curve of the neck, ensuring correct positioning of the patient. This natural neckline alignment allows the neck muscles to relax, preventing tension, neck pain, stiffness, numbness, and muscular inflammation during surgery. Relaxing the neck muscles also helps the shoulder muscles to relax, reducing potential aches and muscular pain during prolonged procedures, which patients frequently complain about post-surgery.

MEMORY FOAM PRONE FACE CUSHION PLACEMENT ILLUSTRATION IN PRONING PROCEDURES



Safer Prone Positioning. Prone Head Supports & MF Cushions

PRONE HEAD SUPPORT MF FACE CUSHIONS

MF PRONE FACE CUSHIONS

Memory Foam Prone Face Cushion


The TIMS MF Prone Face Cushion comes in two sizes, Small and Large, and includes a disposable head measuring tape to assist in choosing the right size.

- Experience a decrease in facial edema and periophthalmitis.
- Lower the risk of periorbital edema, postoperative vision loss, or stroke in your patient.
- Shield your patient's chin from pressure-related injuries.
- Enhance comfort for your patient in the prone position.

references / MF FACE CUSHIONS

Prone Face Cushion Memory Foam	L	1 pc box of 10 pcs	FC-1L-01 FC-1S-10
	S	1 pc box of 10 pcs	FC-1S-01 FC-1S-10

MF PRONE FACE CUSHION SETUP ILLUSTRATION



TIMS

sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Prone Positioning. Prone Head Supports & MF Cushions

PRONE HEAD SUPPORT FACE CUSHIONS


MF SPLIT PRONE FACE CUSHIONS

Memory Foam Split Prone Face Cushion

The **MF Split Prone Face Cushion** allows you to position the patient without without disconnecting their ET Tubing. Pass the tubing through the cut in the side of the cushion and then tape back together.

references / MF CUSHIONS

Memory Foam Split Prone Face Cushion	L	1 pc box of 10 pcs	FC-120-L-SPLIT FC-120-L-SPLIT-10
	S	1 pc box of 10 pcs	FC-120-S-SPLIT FC-120-S-SPLIT-10



NON-DISCONNECTION OF ENDOTRACHEAL TUBING

- *When airway management protocol requires the patient's ET tube to remain connected during safe turning into the prone position, use the Memory Foam SPLIT Face Cushion, which can be placed around the ET tube and then positioned on the patient's face.*
- *Next, slide the St-Airway-Cradle-01 on to the face cushion from the forehead end, around the ET tube. Secure the cradle and Memory Foam Prone Face Cushion with the head strap, and then safely turn the patient into the prone position.*

Safer Prone Positioning. Prone Head Supports & MF Cushions

PRONE HEAD SUPPORT


FACE SUPPORT BASE


Face Support Base

- To achieve greater height between the patient's mouth and the base to accommodate ET Tubing with drainage fittings.
- Disposable prone head support for short surgical procedures.
- Monitor your patient's eyes and tubing during the procedure in the mirror strip.
- Should be used with either the large or small MF Prone Face Cushion.
- **Single use** for infection control

references / SUPPORT

Face Support Base	**[*] sold separately**	FSB-127-02






>> NON-DISCONNECTION OF ET TUBBING WHEN USED WITH THE MEMORY FOAM SPLIT FACE CUSHION



Save time in theatre with fast patient safer prone positioning

TIMS - TECNOLOGIA INOVADORA MAIS SEGURA, LDA



sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Prone Positioning. Torso Support

Safer Prone Positioning. Torso Support


SURGICAL SUPPORTS

TORSO CUSHIONS

Easily **achieve a neutral neckline** and **avoid compression of your patient's abdomen in the prone position**. With no assembly required, the simplicity of the TIMS Torso Support Cushions will reduce the time it takes you to position your patient.


These designs allow for the patient's arms to be positioned forward.

S



CSD-150-S

M LP




CSD-150-M / LOW PROFILE TORSO CUSHION

M XH




CSD-152-M / TORSO CUSHION WITH EXTRA HEIGHT (50 mm)

L



CSD-150-L

L LP




CSD-151-L / LOW PROFILE TORSO CUSHION

Safer Prone Positioning. Torso Support

TORSO SUPPORT

TORSO CUSHIONS

M



CSD-151-M / TORSO CUSHION WITH ARM CHANNELS

Designed to position the patient's arms alongside their body.

/ PRODUCT REFERENCES AND OPTIONS ON THE NEXT PAGES >>

DRAPES & TOPPERS / COVERS

TIMS offers a comprehensive range of accessories designed to extend the usability of single-use surgical supports, while maintaining the highest standards of patient safety and infection control. This assortment includes drapes, memory foam toppers, and wipeable covers. Drapes and memory foam toppers are intended for single-use and are supplied in kits of 10. The wipeable covers, however, are reusable, as they can be effectively sanitised using hospital-grade wipes, ensuring both hygiene and cost-efficiency across multiple procedures.



DRAPE AND MEMORY FOAM TOPPER IN USE

Safer Prone Positioning. Torso Support

TORSO SUPPORT

TORSO CUSHIONS

PRODUCT REFERENCES

Torso Support Memory Foam 44 cm x 34 cm x 14 cm	S	1 pc	CSD-150-S
Low Profile Torso Support Memory Foam 57 cm x 34 cm x 14 cm	M LP	1 pc	CSD-150-M
Torso Support with Extra Height Memory Foam 57 cm x 34 cm x 19 cm	M XH	1 pc	CSD-152-M
Torso Support 59 cm x 54 cm x 23 cm	L	1 pc	CSD-150-L
Low Profile Torso Support 59 cm x 54 cm x 16 cm	L LP	1 pc	CSD-151-L
Torso Support with Arm Channels Memory Foam 66 cm x 53 cm x 18 cm	M	1 pc	CSD-151-M

CAN'T FIND THE PRODUCT YOU'RE LOOKING FOR?

If you need a device that is not currently available, please talk to us. Our designers may be able to collaborate with you and the clinicians to create a new product.

info@tims.pt



Safer Prone Positioning. Torso Support


TORSO SUPPORT DRAPES & TOPPERS

Drape Covers & Memory Foam Toppers

With a combination of drape covers & memory foam toppers you can extend the life of your disposable support. The waterproof drapes will keep the support underneath clean and safe to use again.

- Safely reuse a single use device.
- Improved infection control.
- Offer safe contouring support and improved comfort for the patient by using the memory foam (MF) toppers to help reduce pressure related injury. MF also allows for individual customization of the support to the patient's shape, reducing peak pressure and preventing pressure buildup in certain areas.
- Simple to assemble and easy to use.
- Drapes can be ordered in a kit which includes (x10) Drapes & Memory Foam Toppers, plus the Torso, Dome or Leg Support to be reused.

DRAPE AND MEMORY FOAM TOPPER PLACEMENT ILLUSTRATION



TIMS

sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS


Safer Prone Positioning. Torso Support

TORSO SUPPORT

DRAPES & TOPPERS


references / DRAPE & TOPPER SETS

Memory Foam Torso Support (x10 pcs) Drape & (x10 pcs) Memory Foam Topper (x1 pc) Support & (x10 pcs) Drape & (x10 pcs) Memory Foam Topper	S	CSD-350-S-10 CSD-450-S-03
Memory Foam Low Profile Torso Support (x10 pcs) Drape & (x10 pcs) Memory Foam Topper (x1 pc) Support & (x10 pcs) Drape & (x10 pcs) Memory Foam Topper	M LP	CSD-350-M-10 CSD-450-M-03
Memory Foam Torso Support with Extra Height (x10 pcs) Drape & (x10 pcs) Memory Foam Topper (x1 pc) Support & (x10 pcs) Drape & (x10 pcs) Memory Foam Topper	M XH	CSD-352-M-10 CSD-452-M-03
Torso Support (x10 pcs) Drape & (x10 pcs) Memory Foam Topper (x1 pc) Support & (x10 pcs) Drape & (x10 pcs) Memory Foam Topper	L	CSD-350-L-10 CSD-450-L-03
Low Profile Torso Support (x10 pcs) Drape & (x10 pcs) Memory Foam Topper (x1 pc) Support & (x10 pcs) Drape & (x10 pcs) Memory Foam Topper	L LP	CSD-351-L-10 CSD-451-L-03
Memory Foam Torso Support with Arm Channels (x10 pcs) Drape & (x10 pcs) Memory Foam Topper (x1 pc) Support & (x10 pcs) Drape & (x10 pcs) Memory Foam Topper	M	CSD-351-M-10 CSD-451-M-03





DRAPE & MEMORY FOAM TOPPER IN USE



sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Prone Positioning. Dome Supports

TIMS

sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Prone Positioning. Dome Supports


MF DOME SUPPORTS

MF CHEST & PELVIC

Memory Foam Prone Dome Supports


- The positioning of the Chest and Pelvic Domes to provide improved oxygenation and increased space for the decompression of the abdomen.
- Available in 2 sizes (including kits) for Obese patients the extra large size should be selected.
- Drapes and memory foam toppers available to extend the life of the Dome Supports.

DRAPE AND MEMORY FOAM TOPPER ASSEMBLY ILLUSTRATION



MEMORY FOAM DOME CHEST SUPPORT

The Foam Chest Dome support has been shaped with Memory Foam to provide safe and comfortable support for the female breasts reducing the risk of pressure related injuries and ulcers.




OPTIONS ARE DETAILED IN THE PAGES AHEAD

Safer Prone Positioning. Dome Supports


MF DOME SUPPORTS CHEST

L




DMCC-01L

L LP



DMCC-01LP

XL



DMCC-01XL

/ PRODUCT REFERENCES AND OPTIONS
ON THE NEXT PAGES >>


Safer Prone Positioning. Dome Supports


MF DOME SUPPORTS

CHEST & PELVIC


MEMORY FOAM DOME PELVIC SUPPORT

Foam Pelvic Dome support shaped with Memory Foam to safely position male genitalia, reducing the risk of pressure related injuries and ulcers. Protects male genitalia more safely.





OPTIONS ARE DETAILED IN THE PAGES AHEAD




DRAPE AND MEMORY FOAM TOPPER IN USE

DRAPES & TOPPERS / COVERS


All covering and topper options are available for our leg support cushions, providing the same advantages as our other products. These options ensure that your leg support cushions benefit from enhanced durability, improved fit, and ease of cleaning. By using these covers and toppers, you can maintain a high standard of patient safety and comfort while extending the life of your cushions.

Safer Prone Positioning. Dome Supports


MF DOME SUPPORTS PELVIC



DMPC-01L



DMPC-01LP



DMPC-01XL

/ PRODUCT REFERENCES AND OPTIONS
ON THE NEXT PAGES >>

TIMS


sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS


Safer Prone Positioning. Dome Supports

MF DOME SUPPORTS CHEST & PELVIC

references / DOME CHEST SUPPORT CUSHION

Dome Chest Support Memory Foam 36 cm x 26 cm x 21.5 cm	L	1 pc	DMCC-01L
Low Profile Dome Chest Support Memory Foam 50 cm x 28 cm x 18 cm	L LP	1 pc	DMCC-01LP
Dome Chest Support Memory Foam 60 cm x 30 cm x 24.5 cm	XL	1 pc	DMCC-01XL





DOME CHEST AND PELVIC SUPPORT CUSHIONS

references / DOME PELVIC SUPPORT CUSHION

Dome Pelvic Support Memory Foam 55 cm x 24 cm x 21.4 cm	L	1 pc	DMPC-01L
Low Profile Dome Pelvic Support Memory Foam 60 cm x 27 cm x 19 cm	L LP	1 pc	DMPC-01LP
Dome Pelvic Support Memory Foam 60 cm x 30 cm x 24.5 cm	XL	1 pc	DMPC-01XL

Safer Prone Positioning. Dome Supports

MF DOME SUPPORTS CHEST & PELVIC KITS


references / PRONING KITS

Dome Proning Kit Chest & Pelvic Support kit	L	10 (pcs) DDL-03-DRAPE 10 (pcs) DMCC-TOPPER	DPK-122-01
	L LP	10 (pc) DDL-03-DRAPE 10 (pc) DMPC-TOPPER	DPK-127-01
	XL	10 (pc) DDL-03-DRAPE 10 (pc) DMPC-TOPPER	DPK-126-01





INCREASED SPACE TO ALOW FOR THE DECOMPRESSION OF THE ABDOMEN



CHEST AND PELVIC DOME SUPPORTS SETUP ILLUSTRATION

BETTER POSITIONING

The Foam Chest Dome support has been shaped with Memory Foam to provide safe and comfortable support for the female breasts and the Foam Pelvic Dome support shaped with Memory Foam to safely position male genitalia, reducing the risk of pressure related injuries and ulcers.

Safer Prone Positioning. Dome Supports

MF DOME SUPPORTS

DRAPES & TOPPERS

references / SUPPORT & DRAPE & MF TOPPER SETS

Chest Dome Support & Drape & MF Topper	L	1 (pc) DMCC-01L 10 (pcs) DDL-03-DRAPE 10 (pcs) DMCC-TOPPER	DMCC-450-L-03
	L LP	1 (pc) DMCC-01LP 10 (pcs) DDXL-03-DRAPE 10 (pcs) DMCC-TOPPER	DMCC-450-LP-03
	XL	1 (pc) DMCC-01XL 10 (pcs) DDXL-03-DRAPE 10 (pcs) DMCC-TOPPER-XL	DMCC-450-XL-03

TIMS offers convenient sets that include 10 drapes and 10 memory foam toppers, paired with the Dome Support of your choice. If you already own the support cushion, you can simply select one of our Drape & Topper sets, which provide 10 drapes and 10 toppers, ensuring continued comfort and protection for your patients.

references / DRAPE & MF TOPPER SETS

Dome Chest Support **Drape + MF Topper Set**	L LP	10 (pcs) DDL-03-DRAPE 10 (pcs) DMCC-TOPPER	DMCC-350-L-03
	L	10 (pcs) DDL-03-DRAPE 10 (pcs) DMCC-TOPPER	DMCC-350-LP-03
	XL	10 (pcs) DDXL-03-DRAPE 10 (pcs) DMCC-TOPPER-XL	DMCC-350-XL-03

Safer Prone Positioning. Dome Supports

MF DOME SUPPORTS

WIPEABLE DRAPES & TOPPERS

Reusable/Wipeable Covers

- Extend the lifespan of disposable supports by using our wipeable covers in conjunction with memory foam toppers.
- Replace waterproof single-use drapes with these reusable covers, which help maintain cleanliness and are ready for multiple uses.
- Designed for a better fit, these covers provide enhanced stability and protection.
- Improve infection control and durability while ensuring optimal patient comfort and care.



CLEAN YOUR COVERS WITH HOSPITAL-GRADE WIPES

WIPEABLE / REUSABLE COVERS

With proper sanitisation, wipeable covers can greatly extend the lifespan of surgical supports. Easily cleaned with hospital-grade disinfectants, they allow for safe, repeated use while maintaining high standards of patient care and cost-efficiency and infection control.


Safer Prone Positioning. Dome Supports


MF DOME SUPPORTS

WIPEABLE DRAPES & TOPPERS

references / WIPEABLE COVER

Dome Chest Support **Reusable/Wipeable Cover**	L	1 (pc) wipeable cover	DMCC-RC
	L LP	1 (pc) wipeable cover	DMCC-LP-RC
	XL	1 (pc) wipeable cover	DMCC-XL-RC





REUSABLE / WIPEABLE COVERS IN USE

references / WIPEABLE COVER

Dome Pelvic Support **Reusable/Wipeable Cover**	L	1 (pc) wipeable cover	DMPC-RC
	L LP	1 (pc) wipeable cover	DMPC-LP-RC
	XL	1 (pc) wipeable cover	DMPC-XL-RC



sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Prone Positioning. Arm Supports

Safer Prone Positioning. Arm Supports

ARM SUPPORTS

FOREARM CUSHIONS

Forearm Supports

he Forearm Support safeguards your patient's arms by reducing pressure on the forearms, helping to prevent redness and bruising. It also protects the arms from accidental bumps during procedures.

- The cushion's contoured design allows for safe tilting to position the arms optimally.





references / CUSHIONS

Forearm Cushions (5 pairs) 43 cm	[*] sold in pairs	FASD-150-L

Safer Prone Positioning. Arm Supports

ARM SUPPORTS

FOREARM CUSHIONS





references / CUSHIONS

Forearm Cushions (5 pairs) 56 cm	[*] sold in pairs	FASD-151-L




sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Prone Positioning. Leg Supports

Safer Prone Positioning. Leg Supports


SURGICAL SUPPORTS LEGS

S




LSD-150-S

M



LSD-150-M

L



LSD-150-L


TIMS Leg Support Cushions are available in a range of sizes, designed to accommodate the proning requirements of both adults and children. These cushions can be fitted with custom wipeable reusable covers, enhancing their durability while maintaining high standards of infection control and patient safety. For those who prefer single-use options, we also offer drape and memory foam topper sets.

Safer Prone Positioning. Leg Supports

SURGICAL SUPPORTS LEGS

Adult Leg Support

- Avoid a pressure injury on your patient's knees and allow the patient's feet to sit without pressure on the toes.
- The Prone Leg Support will aid in venous return and prevent pooling of blood.
- Where possible, always ensure that the heels of the feet of the patient are above the surgical site, to prevent air embolisms forming.



LEG SUPPORT CUSHION

The Prone Leg Support ensures proper elevation and hip flexion for your patient, while keeping their feet safely positioned to prevent nerve-related injuries. It supports the shins and knees while allowing the ankles and feet to remain free, reducing the risk of nerve and tissue damage.

OPTIONS ARE DETAILED IN THE PAGES AHEAD

references / LEG SUPPORT CUSHION

Leg Support Cushion 61 cm x 54 cm x 24.5 cm	L	1 pc	LSD-150-L


Safer Prone Positioning. Leg Supports

SURGICAL SUPPORTS

PAEDIATRIC LEG SUPPORTS

Paediatric Leg Support

- The shape of the Prone-Easy Leg Support is designed to promote venous return, preventing blood from pooling in the legs.
- Achieve natural hip and knee flexion when positioning your patient on the Prone-Easy Paediatric Leg Support.
- Allow the ankles and feet to hang free and not come in contact with the table, this can occur when pillows are used due to the pillow collapsing in height over time.
- Where possible, always ensure that the heels of the feet of the patient are above the surgical site, to prevent air embolisms forming.




references / PAEDIATRIC LEG SUPPORT CUSHIONS

Leg Support Cushion 30 cm x 28 cm x 16 cm	S	1 pc	LSD-150-S
Leg Support Cushion 43 cm x 40 cm x 20 cm	M	1 pc	LSD-150-M

Safer Prone Positioning. Leg Supports

ADULT LEG SUPPORTS

DRAPES & TOPPERS



REUSABLE WIPEABLE DRAPE & MEMORY FOAM TOPPER USAGE ILLUSTRATION

references / DRAPE & MF TOPPER SETS

Leg Support **Drape + MF Topper Set**	S	10 (pcs) DCL-03-DRAPE 10 (pcs) LSD-150-S-TOPPER	LSD-350-S-10
	M	10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-M-TOPPER	LSD-350-M-10
	L	10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-L-TOPPER	LSD-350-L-10

references / SUPPORT & DRAPE & MF TOPPER SETS

Leg Support & Drape & MF Topper	S	1 (pc) LSD-150-S 10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-S-TOPPER	LSD-450-S-03
	M	1 (pc) LSD-150-M 10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-M-TOPPER	LSD-450-M-03
	L	1 (pc) LSD-150-L 10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-L-TOPPER	LSD-450-L-03

Safer Prone Positioning. Leg Supports

ADULT LEG SUPPORTS

WIPEABLE DRAPES & TOPPERS


references / WIPEABLE COVERS


Paediatric Leg Support Cushion w/ Resusable Cover	S	1 pc wipeable cover	LSD-150-S-RC
	M	1 pc wipeable cover	LSD-150-M-RC
Leg Support Cushion w/ Resusable Cover **Manufactured in Memory Foam**	L	1 pc wipeable cover	LSD-150-L-RC


references / SUPPORT & WIPEABLE COVER & TOPPER SETS

Leg Support & Wipeable Drape & MF Topper	S	1 (pc) LSD-150-S-RC 10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-S-TOPPER	LSD-550-S-03
	M	1 (pc) LSD-150-M-RC 10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-M-TOPPER	LSD-550-M-03
	L	1 (pc) LSD-150-L-RC 10 (pcs) DCL-03-DRAPE 10 (pair) LSD-150-L-TOPPER	LSD-550-L-03

EXTEND THE LIFESPAN OF THE FOAM LEG SUPPORTS WITH WIPEABLE COVERS







(!) NOTICE: DO NOT AUTOCLAVE. DO NOT SUBMERGE IN WATER



sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Supine Positioning. Supine Head Support


Safer Supine Positioning. Supine Head Support

SURGICAL SUPPORTS BEACH POSITION

Beach Position Head Support

Maintains the patient's head position during surgery.

This device was developed with the help of a Consultant Anaesthetist in Portugal. Suitable for most supine procedures, including shoulder surgery. The head position is easily adjustable with the two angle pads. The memory foam topper provides enhanced customization to the patient's individual head profile and protects the head support and angle pads. The head support can be covered with a drape under the topper for infection control. (Drape not supplied with kit).



BEACH SUPPORT IN USE WITH AN ANGLE PAD FOR BETTER ADJUSTMENT IN A SETTING WERE THE SURFACE IS SET AT AN ANGLE



BEACH SUPPORT USED TO SECURE THE PATIENT'S HEAD IN A FLAT SURFACE

TIMS

sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS


Safer Supine Positioning. Supine Head Support

SURGICAL SUPPORTS BEACH POSITION

references / BEACH HEAD SUPPORT KIT

Beach Position Head Support 190 x 220 x 55 (mm)	**1 box** **7 kits per box**	**1 (pc)** SHSD-01 **1 (pair)** SHSK-AP-02 **20 (pcs)** SHSK-TOPPER	SHSK-20-01

Each Head Support is supplied with a (x1) Memory Foam Topper and 20 Angle Pads for improved adjustability.



ADJUST THE PATIENT'S HEAD IN THE SUPINE POSITION WITH TIMS'S BEACH POSITION HEAD SUPPORT



sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

Safer Patient Care. Accessories





Safer Patient Positioning. Accessories

ACCESSORIES PROTECTIVE STRAPS

Protective Straps

- Available in 6 different lengths and widths.
- Made with brushed nylon and is soft and kind to the skin.
- Strong and reliable.
- Easily cut to any size.
- Full adjustable VELCRO® hook & loop locking system.
- Single use for superior infection control.








TIMS NYLON STRAPS WITH VELCRO® PROVIDE A SAFE AND RELIABLE SOLUTION FOR MEDICAL STAFF TO SECURELY POSITION PATIENTS DURING PROCEDURES

Safer Patient Positioning. Accessories


ACCESSORIES PROTECTIVE STRAPS



Arm Straps
SCS-01



Body Straps
SCS-02



Knee Straps
SCS-03

references /BRUSHED NYLON PROTECTIVE STRAPS

Arm Strap 800 mm x 40 mm	**1 Velcro® fastener**	**50 un** per box	SCS-01
Body Strap 1600 mm x 100 mm	**1 Velcro® fastener**	**14 un** per box	SCS-02
Body Strap 1600 mm x 100 mm	**2 Velcro® fasteners**	**14 un** per box	SCS-06
Body Strap 900 mm x 75 mm	**2 Velcro® fasteners**	**14 un** per box	SCS-07
Knee Strap 1500 mm x 75 mm	**1 Velcro® fastener**	**14 un** per box	SCS-03
Shoulder Stop Strap 1000 mm x 75 mm	**2 Velcro® fasteners**		SCS-09

Safer Patient Positioning. Accessories

SAFER PATIENT POSITIONING — MIRRORS

Mirrors



Mirror		
260 x 110 (mm)	1 un	PLM-01
180 x 35 (mm)	1 un	M-180X35
210x130 (mm)	1 un	M-210X130

SAFER PATIENT POSITIONING — HELMET ANGLE PADS

Helmet Angle Pads

Adjust the height and angle of a prone head support to achieve the optimum natural neutral neckline for the patient.



Helmet Angle Pads [*] 1 kit	1 (pair) Angle Pad Disposable 1 (pc) Flat Pad Disposable	HAP-01D

Safer Eye Protection. Accessories

SAFER CLINICIAN EYE PROTECTION — EYE SHIELDS

Eye Shield Frames & Lenses

Our Eye Shield Frames & Lenses have proven extremely popular for protection against the splash and splatter of blood borne pathogens and other potentially infectious materials.



Eye Shield Lenses 250 (pcs) per pack	9210-250
Eye Frames 100 (pcs) per pack	9211-100

SAFER PATIENT EYE PROTECTION — EYE PROTECTION

Eye Protection

Opti-Gard Patient Eye Protection provides non-allergenic, self-adhesive sponge cushions that reduce the evaporation of liquid from the surface of the eye, which can happen in absence of spontaneous blinking, causing the cornea to dry out.



Opti-Gard 30 (pcs) per box	OPTI-30A

Safer Patient Health Protection. Accessories

SAFER PATIENT HEALTH PROTECTION

PROTECTIVE MITTS

Protective Mitts

The comfortable and ergonomic design aids in reducing emotional and physical discomfort which an agitated patient may feel, when restrained.



Protective Mitts MITTS-02-L

[*] 5 pairs per box

NEED FURTHER ASSISTANCE?

SCAN ME



If you have any questions or need assistance after reviewing our catalogue, please don't hesitate to reach out. We're here to help with any inquiries you may have at **info@tims.pt**.

TIMS

SUMEDEX S.A.
C/TORRETA DE MIRAMAR 12 BAJO,
46020 VALENCIA

sumedex, s.a.
SUMINISTROS MEDICOS EXCLUSIVOS

E-mail: sumedex@sumedex.es

SCAN ME